

echoes Depot





Medstrat has created **echoes Depot**, an entirely automated, off-site image backup and archiving solution.

echoes Depot

Just like physical film is susceptible to disasters, digital images on servers and hard drives can also be lost if equipment fails. That's why every Medstrat PACS system provides the ability to backup images to another server, external hard drives or DVD media. As your image library grows, your backup and archive storage needs will also grow which means you're always going to need more and more data storage.

In response to customer requests, Medstrat created **echoes Depot**, an entirely automated, off-site image backup and archiving solution. Now customers no longer have to worry about backups and archiving and dealing with additional storage media or hardware. The **echoes Depot** solution takes care of it all. Think of **echoes Depot** as a really huge **echoes Server** that's located at the Medstrat datacenter. The system does all the work for you in the background so you don't have to worry about safeguarding your data.



The **echoes** Depot solution allows unlimited, scalable storage and is the best disaster recovery protection for your PACS system available in the market today.



Making Disaster Recovery Far Less Disastrous

When the **echoes** Server in your office is connected remotely to **echoes** Depot, all of the patient images and information on your **echoes** Server are backed up and archived in real-time 24/7. The **echoes** Depot solution allows unlimited, scalable storage and is the best disaster recovery protection for your PACS system available in the market today. The **echoes** Depot will literally contain a mirror image of your **echoes** Server along with all of your backed up and archived images so that if the **echoes** Server in your office was ever to fail, you can easily log into **echoes** Depot remotely and access all of your images just as if your **echoes** Server was still up and running.



MEDSTRAT

Medstrat is Orthopedic PACS

1901 Butterfield Road • Suite 600 • Downers Grove, IL 60515
medstrat.com